



SPECTROLAB



Newbury, Berkshire. RG14 1JB UK. Offices in UK, EU, Dubai

sales@spectrolab.eu www.spectrolab.co.uk

News Release NL 338
September 2019



Platinum + other elements Crucibles and Labware



-4#

Spectrolab are now supplying its customers with Platinum tailor-made products for every need. They include Platinum crucibles and accessories for XRF analysis, microanalysis, burning and evaporation (whenever resistance to high temperatures, chemical attack by acids and fusion mixtures are required, The range includes mini-crucibles, moulds (with or without pouring lips, pins on a supporting ring, reinforced rims, casting moulds, boats, tongs with shot- blasted platinum shoes, mixers, triangles, and electrodes. All are manufactured in pure platinum or alloys with Gold, Rhodium, Iridium, Palladium and other precious metals.



SPECTROLAB



Newbury, Berkshire. RG14 1JB UK. Offices in UK, EU, Dubai
sales@spectrolab.eu www.spectrolab.co.uk