



SPECTROLAB



Newbury, UK. Offices in UK, EU, Dubai
sales@spectrolab.eu sales@spectrolab.co.uk www.spectrolab.co.uk

News Letter 440

June 2022

Best Value... Ever

Prices from just £1850 Why pay more ?

“ Multi element flame photometers”

Na, K, Ca, Li, Ba, Rb, Cs, Sr

**For industry, hospitals, environment control, colleges,
universities, water analysis... anyone**

Sensitivity 0.1ppm



Model S20-8 8 Element FP

The Spectrolab S20-8 multi element flame photometer is an extremely accurate PC or microprocessor-controlled instrument designed for both environmental and medical applications.

Features

- 8 element simultaneous multichannel analysis
- 4 Models available, 8 element, 5 element 3 element 2 element
- Simultaneous analysis with a single aspiration
- Stores multiple calibrations
- Auto ignition, auto gas shut down with power failure
- Auto correction for interferences
- Microprocessor data storage
- Connect a computer for software control
- Built in printer or connect a printer
- USB and RS232 interface

Delivery One week

Competitors ? With this price and performance.....NONE

The S20 - 8 *The only game in town*